

MB90 and AQUALAB 3 Setup Checklist

Oct 1, 2026 · @Arnaud

Print this sheet and tick each box as you go. If a step fails, run `1_check_setup.py` from the kit's Troubleshooting folder and note what it says in the Notes section at the bottom.

MB90 moisture analyzer

Cable

- ☐ USB-A to USB-B cable: flat end in the laptop, square end in the MB90

Settings on the MB90 (wrench icon → Device settings)

- ☐ Peripheral settings: Device Port = **COM4**, 9600 baud, 8 data bits, parity None, 1 stop bit, handshake None
- ☐ Print settings: **Print measurement result automatically = ON** and **Print adjustment result automatically = ON**
- ☐ Date-Time settings: date and time are correct

Connect and test (close the logger first)

- ☐ Run `2_find_ports.py`, then plug the cable into the **square USB-B port** on the back of the MB90
- ☐ New port noted (ignore COM4 and COM5, the laptop's Bluetooth ports): COM _____
- ☐ Run `3_test_port.py COM__`, choose **1. Listen**, press **Print** on the MB90: verdict says **LOOKS GOOD**
- ☐ `lab_logger_config.json`: MB90 "port" set to the new port
- ☐ Start the logger, type `status`: **MB90: connected**
- ☐ Run a short test: logger shows **MB90: saved**, and `data\MB90\MB90_2026-10.csv` matches the MB90 screen

If you see	Do this
No new port	Try another laptop USB port; check Device Manager → Other devices for a yellow icon (needs OHAUS driver); try another cable

If you see	Do this
Garbled text	In <code>3_test_port.py</code> choose 5. Find the right baud rate
Nothing when Print is pressed	Recheck Device Port = COM4 on the MB90; unplug and replug

AQUALAB 3 / SKALA Control

The AQUALAB 3 doesn't plug into the laptop. For now, readings are downloaded from SKALA Control into a folder the logger watches. Downloading needs internet.

Update the logger (one time)

- ☐ New `lab_logger.py` copied over the old one in `C:\LabLogger`
- ☐ **AQUALAB3** entry added to `lab_logger_config.json` (connection guide, Step 4A)
- ☐ `1_check_setup.py` shows `AQUALAB3: watch folder` and no JSON errors

Test a download

- ☐ Signed in to SKALA Control while online
- ☐ One day's readings downloaded as CSV or .xlsx
- ☐ File saved into `C:\LabLogger\data\AquaLab Exports`
- ☐ Logger shows **AQUALAB3: imported**
- ☐ `data\AQUALAB3\AQUALAB3_2026-10.csv` matches the readings in SKALA
- ☐ Same day downloaded again: logger says **skipped** and adds no new rows

Use SKALA's own **batch** column for AQUALAB 3 lot numbers. The logger's `lot` column only shows what was typed when the file was imported.

To collect for automatic SKALA fetching

Send these to Claude so the logger can pull AQUALAB 3 readings by itself whenever the laptop is online. All of them come from Chrome's DevTools: press **F12** on skalacontrol.io and open the **Network** tab.

- ☐ **The request:** filter for `get_readings_by_test`, right-click it → Copy → **Copy as cURL (bash)**, and paste it into a text file

- ☐ **Secrets removed:** in that file, replace the values after **Authorization:** and **Cookie:** , and any token, session or password fields, with **XXXX**
- ☐ **Sample response:** from the request's **Response** tab, one or two readings copied out
- ☐ Optional: **the login request** (sign out, clear the Network log, sign in, and copy the request named like **login** or **token** , with secrets removed the same way)

Never send working passwords or tokens in chat. The logger will read them from the config file on the laptop.

Final full test

- ☐ Logger running; typed **op <name>** and **lot TEST**
- ☐ MB90: one test logged, CSV matches the screen
- ☐ pH meter: one reading printed, CSV matches the screen
- ☐ LS177C: one scan exported, CSV matches the software
- ☐ AQUALAB 3: one download imported, CSV matches SKALA Control

Notes

Date	Instrument	What happened / port / error